



Transmembrane Protein MDS

- 1.) Open .pdb file inside of YASARA (.sce & .job files will not work for this purpose)
- 2.) Select the target by clicking Options-> Macro & Movie-> Set Target (this is when you will select the model that you currently have open)




- 3.) Now tell YASARA which macro to run by clicking Options-> Macro & Movie-> Play Macro-> md_runmembrane.mcr




Transmembrane Protein MDS

- 4.) YASARA will now scan your protein for exposed hydrophobic regions that suggest potential transmembrane regions




- 5.) A suggested membrane will be placed on your protein.



- a. You can edit location and orientation of your membrane by clicking shift twice to control Object 3

Transmembrane Protein MDS

- b. Select Continue when content with membrane location and orientation



- c. IMPORTANT: Yasara is not always accurate. Be sure to compare predicted membrane location with outside sources.

- 6.) YASARA will now build the actual membrane using default conditions. These conditions can be edited in the actual macro document. (current options are phosphatidyl-ethanolamine, -choline, and -serine)
- 7.) YASARA will slowly build out the membrane inside of the simulation cell
- 8.) Occasionally there will be an error stating that the membrane exploded
 - a. This usually means that your membrane location was not optimal. Trial and error is the best way to fix this issue.
- 9.) Once the membrane has been completely built a 250 picosecond molecular dynamic simulation will occur to stabilize your protein in the membrane environment.
- 10.) You should now have a protein embedded in a lipid transmembrane environment.